


1. Introduction


Neuroinflammation is increasingly linked to the etiology and progression of Alzheimer's disease (AD). Microglia, the critical stewards of the brain's immune microenvironment, display phenotypic and state diversity in AD as shown in single-cell transcriptomic and proteomic analyses on tissue homogenates¹. These approaches lack spatial information, which is crucial to interpreting microglial involvement in neuroinflammation. We thus set out to develop a workflow that provides both, deep biomarker readouts and excellent spatial resolution. To this end, we developed a **PhenoCycler®-Fusion (PCF)** workflow that is compatible with formalin-fixed human brain tissue. The PCF is a novel spatial biology platform that generates ultrahigh-plex spatial data for millions of cells on a single tissue section within a matter of hours. Using this technology, we characterized the immune microenvironment – with emphasis on microglia – in normal and AD brain tissues. Our findings suggest an enrichment of certain microglial subtypes in AD brain tissues and in proximity to Aβ plaques. In its sum, this study provides (1) a method for ultrahigh-plex immunofluorescent analyses of AD pathology and new insights into the inflammatory microenvironment of AD.

2. PhenoCycler-Fusion Workflow




The PhenoCycler-Fusion workflow consists of iterative cycles of labelling, imaging and removing fluorescent reporters. In each imaging cycle, three fluorescent reporters are assayed to their corresponding barcode-conjugated antibodies and imaged via standard fluorescent optics. Thereafter, the three reporters are removed, and a new cycle images additional reporters. The process is fully automated, and data are acquired **across whole slides at single-cell resolution**.


3. Comprehensive Mapping of Alzheimer's Pathology in Human Brain Tissue



4. Stardist Nuclear Segmentation for Neuronal Cell Classification




(A) Fluorescent image showing major cell-type defining markers. (B) Spatial map showing cell classification results for 110,088 DAPI-based segmentation masks, obtained using Stardist. GM: Grey Matter, WM: White Matter (WM).




The Stardist cell masks failed to capture cellular processes in astrocytes and microglia (see GFP & IBA1 on left). To correctly segment individual microglia, we thus developed **OTSM** (pronounced 'Awesome'), which stands for Otsu's Thresholding based Segmentation and Merge. This method obtains a binary mask around microglial cells using Otsu's thresholding (2), followed by a merging step to merge disconnected processes using a predetermined distance cutoff. **Box 5** summarizes our results.

5. OTSM ('Awesome') Segmentation Accurately Registers Microglia



6. Inflammatory State of Microglia Depends on Proximity to Aβ Plaques



7. Conclusions and Outlook

- More than five million Americans are living with Alzheimer's Disease, but our understanding of the etiology and progression of this disease remains limited.
- We developed a comprehensive Spatial Biology workflow aimed at uncovering the inflammatory biology of AD in human patient samples *in situ*.
- Our work encompasses the development of a custom antibody panel, an imaging workflow, as well as a novel bioinformatic analysis method.
- Deployment of this workflow on post-mortem AD tissues allowed us to study different microglial cell populations, according to biomarker profiles and spatial distribution; in doing so, we identified different area and disease-associated microglial subpopulations in the cerebral cortex.
- Our workflow will allow large-scale unsupervised analyses of the AD inflammatory microenvironment that is needed to better characterize neurodegenerative disease.

1. Wang H (2021). Microglia Heterogeneity in Alzheimer's Disease: Insights From Single-Cell Technologies. *Front. Synaptic Neurosci.*
2. N. Otsu (1979). A Threshold Selection Method from Gray-Level Histograms, *IEEE Trans. on Sys, Man, & Cybernetics*

This study was supported by Alzheimer's Association Grant #AARGD-22-92829 to BA.

